

Product datasheet for **MP205966**

Hsp90b1 Mouse qPCR Primer Pair (NM_011631)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22027
Forward Sequence:	GTTTCCCGTGAGACTCTTCAGC
Reverse Sequence:	ATTCGTGCCGAATCCTTCAG
ACCN:	<u>BC010445</u> , <u>BC011439</u> , <u>NM_011631</u> , <u>NM_011631.1</u>
UniProt ID:	<u>P08113</u>
Synonyms:	endoplasmin; ERp99; gp96; GRP94; TA-3; Targ2; Tra-1; Tra1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)