

## Product datasheet for **MP205967**

### Hsp90ab1 Mouse qPCR Primer Pair (NM\_008302)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15516
Forward Sequence:	CCTGAAGGTCATCCGCAAGAAC
Reverse Sequence:	GGCGTCGGTTAGTGGAATCTTC
ACCN:	<a href="#">BC088985</a> , <a href="#">NM_008302</a> , <a href="#">NM_008302.1</a> , <a href="#">NM_008302.2</a> , <a href="#">NM_008302.3</a> , <a href="#">BC002167</a> , <a href="#">BC044888</a> , <a href="#">BC049951</a> , <a href="#">BC054457</a>
UniProt ID:	<a href="#">P11499</a>
Synonyms:	90kDa; AL022974; C81438; Hsp84; Hsp84-1; Hsp90; Hspcb
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)