

Product datasheet for **MP205977**

Herc4 Mouse qPCR Primer Pair (NM_026101)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67345
Forward Sequence:	CCAGGCAAAACTACTCTGCTCC
Reverse Sequence:	GGATTCACAGACTCGATCACTGG
ACCN:	BC043082 , BC060033 , NM_026101 , NM_026101.1 , NM_026101.2 , NM_026101.3 , NM_026101.4 , BC026855 , BC047157 , BM249989 , BQ899503
UniProt ID:	Q6PAV2
Synonyms:	1700056O17Rik; 4921531D01Rik; 9530080M15Rik; mKIAA1593
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)