

Product datasheet for **MP205987**

HegI Mouse qPCR Primer Pair (NM_175256)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	77446
Forward Sequence:	TGCTCAGAGCACCTCACAGTCA
Reverse Sequence:	GCCAGTGAGGTCCACAATGCAT
ACCN:	<u>BC076596</u> , <u>NM_175256</u> , <u>NM_175256.1</u> , <u>NM_175256.2</u> , <u>NM_175256.3</u> , <u>NM_175256.4</u> , <u>NM_175256.5</u> , <u>BB639887</u> , <u>BC038614</u> , <u>BC072651</u> , <u>BE862445</u>
Synonyms:	4632417D23Rik; 5530401I02Rik; 9530025L16Rik; BC038614; Gm629; OTTMUSG00000027315
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

