

Product datasheet for **MP205994**

Jpt2 Mouse qPCR Primer Pair (NM_198937)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52009
Forward Sequence:	CCTCTGGCACACCCAAACAAGC
Reverse Sequence:	GCTTCCTTTGTCACCTCTGCTCTC
ACCN:	<u>BC057024</u> , <u>BC080740</u> , <u>NM_198937</u> , <u>NM_198937.1</u> , <u>NM_198937.2</u> , <u>BC055351</u> , <u>BK001574</u>
UniProt ID:	<u>Q6PGH2</u>
Synonyms:	2810430B18Rik; AL022629; C16orf34; D17Ertd441e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)