

Product datasheet for **MP206024**

Hs3st3b1 Mouse qPCR Primer Pair (NM_018805)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54710
Forward Sequence:	TTTCGACCGCAGCTACCACAAG
Reverse Sequence:	AAGTAGCTGGGCGTCTTCTCCA
ACCN:	<u>BC116733</u> , <u>BC116743</u> , <u>NM_018805</u> , <u>NM_018805.1</u> , <u>NM_018805.2</u>
UniProt ID:	<u>Q9QZS6</u>
Synonyms:	3-OST-3B; 3Ost3b; 3OST3B1; AU041882; AW536289; Hs3st3b; m3-OST-3B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)