

## Product datasheet for **MP206037**

### Hepacam Mouse qPCR Primer Pair (NM\_175189)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72927
Forward Sequence:	GTCTACCTGCTTTCATCCAGC
Reverse Sequence:	TTGTCGCTGCTGGTGCTACTGT
ACCN:	<u><a href="#">NM_175189</a></u> , <u><a href="#">NM_175189.1</a></u> , <u><a href="#">NM_175189.2</a></u> , <u><a href="#">NM_175189.3</a></u> , <u><a href="#">NM_175189.4</a></u> , <u><a href="#">BC139057</a></u> , <u><a href="#">BC082537</a></u> , <u><a href="#">BC139058</a></u> , <u><a href="#">BY753058</a></u>
UniProt ID:	<u><a href="#">Q640R3</a></u>
Synonyms:	2900042E01Rik; Glialcam; Hepn1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)