

Product datasheet for **MP206061**

Hnrnpc Mouse qPCR Primer Pair (NM_016884)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15381
Forward Sequence:	GCTCTGTGCATAAAGGCTTTGCC
Reverse Sequence:	CTCGGTTCACTTTTGGCTCTGC
ACCN:	<u>BC004706</u> , <u>BC095922</u> , <u>NM_016884</u> , <u>NM_016884.1</u> , <u>NM_016884.2</u> , <u>NM_016884.3</u> , <u>BB283438</u> , <u>NM_016884.4</u>
UniProt ID:	<u>Q9Z204</u>
Synonyms:	AL022939; D14Wsu171e; hnrnp-C; hnRNPC1; hnRNPC2; Hnrpc; Hnrpc1; Hnrpc2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)