

Product datasheet for **MP206063**

Hnrnpdl Mouse qPCR Primer Pair (NM_016690)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50926
Forward Sequence:	ACATACACAGACGAAGAGCCAGT
Reverse Sequence:	GCCTCCTTTCTGTTGTTGCTGC
ACCN:	<u>BC021374</u> , <u>NM_016690</u> , <u>NM_016690.1</u> , <u>NM_016690.2</u> , <u>NM_016690.3</u> , <u>NM_016690.4</u> , <u>BM950445</u> , <u>BQ174478</u>
Synonyms:	AA407431; AA959857; D5Ert650e; D5Wsu145e; hnRNP-DL; hnRNP DL; Hnrnpdl; Hnrpd; JKTBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

