

Product datasheet for **MP206067**

Hnrnp3 Mouse qPCR Primer Pair (NM_001079824)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	432467
Forward Sequence:	GAGCTGGTGATGCAAGTTCAGG
Reverse Sequence:	TTGCTCTGCCATCGGCTCCAAT
ACCN:	NM_001079824 , NM_001079824.1 , BC049652 , BC147496 , BC160206 , BQ443117 , BU946493 , BX514270
Synonyms:	AA693301; A1666703; Hnrp3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

