

Product datasheet for **MP206093**

Hdlbp Mouse qPCR Primer Pair (NM_133808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110611
Forward Sequence:	AGCGACACTGTTGTCATCAGGG
Reverse Sequence:	GGTATTCTGGCTTGCGTCGGAT
ACCN:	<u>BC021765</u> , <u>BC023806</u> , <u>BC025648</u> , <u>BC027779</u> , <u>BC027788</u> , <u>BC035301</u> , <u>NM_133808</u> , <u>NM_133808.1</u> , <u>NM_133808.2</u> , <u>NM_133808.3</u> , <u>NM_133808.4</u> , <u>NM_133808.5</u> , <u>BB854571</u> , <u>BC008524</u> , <u>BC027067</u> , <u>BC029338</u> , <u>BC046545</u>
UniProt ID:	<u>Q8VDJ3</u>
Synonyms:	1110005P14Rik; AA960365; A1118566; D1ErtD101e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

