

## Product datasheet for **MP206096**

### Hmga1 Mouse qPCR Primer Pair (NM\_016660)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15361
Forward Sequence:	AGTGCCAACTCCGAAGAGACCT
Reverse Sequence:	TTGGTTTCCTCCCTGGAGCTGT
ACCN:	<a href="#">BC008125</a> , <a href="#">BC013455</a> , <a href="#">BC093496</a> , <a href="#">NM_016660</a> , <a href="#">NM_016660.1</a> , <a href="#">NM_016660.2</a> , <a href="#">NM_016660.3</a> , <a href="#">BM937661</a> , <a href="#">BY165675</a> , <a href="#">BY270784</a>
UniProt ID:	<a href="#">P17095</a>
Synonyms:	AL023995; Hmg; HMG-I(Y; Hmga; Hmga1a; Hmga1b; Hmgi; HMGI(Y; Hmgii; Hmgy
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)