

Product datasheet for **MP206099**

Hmgb2 Mouse qPCR Primer Pair (NM_008252)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	97165
Forward Sequence:	GATGTGGTCTGAGCAATCTGCC
Reverse Sequence:	CCTGCTTCACTTTTGCCCTTG
ACCN:	BC002050 , BC046759 , BC083108 , NM_008252 , NM_008252.1 , NM_008252.2 , NM_008252.3 , BX516410 , BY009449
UniProt ID:	P30681
Synonyms:	C80539; HMG-2; Hmg2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)