

Product datasheet for **MP206116**

Hrh2 Mouse qPCR Primer Pair (NM_001010973)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15466
Forward Sequence:	TCTCTGCTGTCCAGGAGCCAAA
Reverse Sequence:	TCAGGAGACAACCTTGAGTGACTC
ACCN:	NM_001010973 , NM_001010973.1 , NM_001010973.2 , BB647802 , BC101958 , BC103552 , BC103553 , BC103554 , BC165958
UniProt ID:	P97292
Synonyms:	H2R; HH2R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)