

Product datasheet for **MP206134**

H2-Ab1 Mouse qPCR Primer Pair (NM_207105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14961
Forward Sequence:	GTGTGCAGACAACTACGAGG
Reverse Sequence:	CTGTCACTGAGCAGACCAGAGT
ACCN:	BC008168 , BC010322 , BC057998 , NM_207105 , NM_207105.1 , NM_207105.2 , NM_207105.3 , BB056841 , BY205244
UniProt ID:	P14483
Synonyms:	Abeta; A1845868; H-2Ab; H2-Ab; I-Abeta; Ia-2; Ia2; IAb; Rmcs1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

