

Product datasheet for **MP206145**

H2-MI0.1 Mouse qPCR Primer Pair (NM_013544)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	14985
Forward Sequence:	AAGGCTTCTTAGATACCTGGACC
Reverse Sequence:	GGGTGACATTCCTTCAGGTCT
ACCN:	<u>BC119458</u> , <u>NM_013544</u> , <u>NM_013544.1</u> , <u>NM_013544.2</u> , <u>NM_013544.3</u> , <u>BC104320</u> , <u>BC104321</u> , <u>BC152395</u>
Synonyms:	9.5H
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

