

Product datasheet for **MP206171**

Selenos Mouse qPCR Primer Pair (NM_024439)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109815
Forward Sequence:	CATTGTCATCCAGAGGCTCTCC
Reverse Sequence:	CTGCTAAAGCCTCTTGCCGCTT
ACCN:	<u>BC011091</u> , <u>NM_024439</u> , <u>NM_024439.1</u> , <u>NM_024439.2</u> , <u>NM_024439.3</u> , <u>BC096566</u> , <u>BY335100</u>
UniProt ID:	<u>Q9BCZ4</u>
Synonyms:	1500011E07Rik; C78786; H-4; H-47; H4; H47; Se; Sels; Seps1; V; Vimp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

