

## Product datasheet for **MP206177**

### Hlf2 Mouse qPCR Primer Pair (NM\_015786)

#### Product data:

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	50708
<b>Forward Sequence:</b>	TGGTGCTGCCACCCAAAAAG
<b>Reverse Sequence:</b>	GCTTTAGACGCGCTCTTGACCT
<b>ACCN:</b>	<u><a href="#">BC089604</a></u> , <u><a href="#">NM_015786</a></u> , <u><a href="#">NM_015786.1</a></u> , <u><a href="#">NM_015786.2</a></u>
<b>UniProt ID:</b>	<u><a href="#">P15864</a></u>
<b>Synonyms:</b>	0610008C09Rik; H1-2; H1.2; H1c; H1v; H1var1; His1a; Hist1h; Hist1h1c
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.

