

Product datasheet for **MP206180**

Hlf6 Mouse qPCR Primer Pair (NM_010377)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107970
Forward Sequence:	GCTCAGTAAGAAGGCGGCTTCT
Reverse Sequence:	GGCTTCTTGACAGCCTTGGTCT
ACCN:	BC119305 , BC119307 , NM_010377 , NM_010377.1 , NM_010377.2 , NM_010377.3 , BV210386 , BX630865
UniProt ID:	Q07133
Synonyms:	H1-6; H1.6; H1ft; H1t; Hist1h; Hist1h1t
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

