

Product datasheet for **MP206215**

H4c1 Mouse qPCR Primer Pair (NM_178192)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	326619
Forward Sequence:	AACATCCAGGGCATCACCAAGC
Reverse Sequence:	GTTCTCCAGGAACACCTTCAGC
ACCN:	<u>BC119241</u> , <u>BC119243</u> , <u>NM_178192</u> , <u>NM_178192.1</u> , <u>NM_178192.2</u>
UniProt ID:	<u>P62806</u>
Synonyms:	H4c2; H4c3; H4c4; H4c6; H4c8; H4c9; H4c11; H4c12; H4c14; H4f16; Hist1; Hist1h4a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)