

Product datasheet for **MP206225**

H2acI8 Mouse qPCR Primer Pair (NM_013549)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15267
Forward Sequence:	GCGACAACAAGAAGACGCGCAT
Reverse Sequence:	ATCGTCACTTTGCCAGCAGCT
ACCN:	<u>BC010564</u> , <u>BC062255</u> , <u>BC089519</u> , <u>NM_013549</u> , <u>NM_013549.1</u> , <u>NM_013549.2</u> , <u>BC010564.1</u> , <u>BC053966</u> , <u>BC064002</u> , <u>BC147778</u> , <u>BC147781</u> , <u>BC160305</u>
UniProt ID:	<u>Q6GSS7</u>
Synonyms:	AI448581; H2a-61; H2a-614; H2a-615; His; Hist2; Hist2h; Hist2h2aa1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

