

Product datasheet for **MP206231**

H3c13 Mouse qPCR Primer Pair (NM_178215)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319154
Forward Sequence:	CTGATCCGCAAGCTGCCGTTC
Reverse Sequence:	CACAGTTGGTGTCTCGAACA
ACCN:	<u>NM_178215</u> , <u>NM_178215.1</u> , <u>NM_178215.2</u> , <u>BC139251</u> , <u>BC139252</u> , <u>BC156107</u> , <u>BC156965</u>
UniProt ID:	<u>P84228</u>
Synonyms:	H3-615; H3-616; H3c2; H3c3; H3c4; H3c6; H3c7; H3c14; H3c15; Hist2h; Hist2h3b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)