

Product datasheet for **MP206256**

Kat5 Mouse qPCR Primer Pair (NM_178637)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81601
Forward Sequence:	TGAGCGTGAAGGACATCAGTGG
Reverse Sequence:	TTAAGTCCAGCCGCTCGTGAGT
ACCN:	BC110675 , BC129968 , NM_178637 , NM_178637.1 , NM_178637.2 , BC022768 , BC129967 , BY269987
UniProt ID:	Q8CHK4
Synonyms:	AI839539; CPLA2; Htatip; Htatip1; PLIP; Tip55; Tip60
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)