

## Product datasheet for **MP206262**

### Hnf1a Mouse qPCR Primer Pair (NM\_009327)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	21405
Forward Sequence:	AGAGACCTTGGTGGAGGAGTGT
Reverse Sequence:	GGCAAACCAAGTTGTAGACACGC
ACCN:	<u><a href="#">BC080698</a></u> , <u><a href="#">NM_009327</a></u> , <u><a href="#">NM_009327.1</a></u> , <u><a href="#">NM_009327.2</a></u> , <u><a href="#">NM_009327.3</a></u> , <u><a href="#">BC013454</a></u>
UniProt ID:	<u><a href="#">P22361</a></u>
Synonyms:	AI323641; Hnf-1; HNF1; HNF1-alpha; Hnf1 alpha; HNF1[a]; LFB1; Tcf-1; Tcf1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)