

Product datasheet for **MP206307**

Hmbox1 Mouse qPCR Primer Pair (NM_177338)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219150
Forward Sequence:	GCCTCAGAAGAGGACCTAGATG
Reverse Sequence:	TCCTCCGATTGGCAAGAAAGGC
ACCN:	<u>BC051457</u> , <u>NM_177338</u> , <u>NM_177338.1</u> , <u>NM_177338.2</u> , <u>NM_177338.3</u> , <u>NM_177338.4</u> , <u>NM_177338.5</u> , <u>NM_177338.6</u> , <u>BC002212</u>
UniProt ID:	<u>Q8BJA3</u>
Synonyms:	AI451877; AI604847; F830020C16Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)