

Product datasheet for **MP206313**

Hipk3 Mouse qPCR Primer Pair (NM_010434)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15259
Forward Sequence:	AGATCCGACCTGGAGTCCTTTC
Reverse Sequence:	GCCTTCAGAAGTTAGTGTTGCGG
ACCN:	<u>BC131664</u> , <u>NM_010434</u> , <u>NM_010434.1</u> , <u>NM_010434.2</u> , <u>BC014856</u> , <u>BC050118</u> , <u>BC131663</u>
UniProt ID:	<u>Q9ERH7</u>
Synonyms:	DYRK6; FIST3; mir-1902; PKY
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)