

## Product datasheet for **MP206346**

### Hip1 Mouse qPCR Primer Pair (NM\_146001)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	215114
Forward Sequence:	ATGCCCGAGTTGGAGAAGGAAC
Reverse Sequence:	CTTCACCTGCTGGCACATGGAT
ACCN:	<u><a href="#">BC017516</a></u> , <u><a href="#">NM_146001</a></u> , <u><a href="#">NM_146001.1</a></u> , <u><a href="#">NM_146001.2</a></u> , <u><a href="#">BY309818</a></u>
UniProt ID:	<u><a href="#">Q8VD75</a></u>
Synonyms:	2610109B09Rik; A930014B11Rik; E130315I21Rik; mKIAA4113
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)