

Product datasheet for **MP206352**

Hapln1 Mouse qPCR Primer Pair (NM_013500)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	12950
Forward Sequence:	TACTGCTGTGGTGGCATTGGAG
Reverse Sequence:	GGTCAAAGGAAGCAATCACTGCG
ACCN:	<u>BC066853, NM_013500, NM_013500.1, NM_013500.2, NM_013500.3, NM_013500.4</u>
UniProt ID:	<u>Q9QUP5</u>
Synonyms:	BB099155; CLP; Crtl1; Crtl1l; LP; LP-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)