

Product datasheet for **MP206353**

Hapln2 Mouse qPCR Primer Pair (NM_022031)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73940
Forward Sequence:	CAATGGCATTGAGGACGAGAGC
Reverse Sequence:	CGAAGTAATTGAACTGGTAGCGG
ACCN:	NM_022031 , NM_022031.1 , NM_022031.2 , BC138584 , BC023911 , BC138585 , BE651789 , BY122070
UniProt ID:	Q9ESM3
Synonyms:	4930401E20Rik; Brall
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

