

Product datasheet for **MP206376**

Haghl Mouse qPCR Primer Pair (NM_026897)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68977
Forward Sequence:	CACATTCTGCCAGGATTGGCTG
Reverse Sequence:	CTGGCGTCAAGAGGCATCGCA
ACCN:	BC030466 , BC083322 , NM_026897 , NM_026897.1 , NM_026897.2 , NM_026897.3 , BC005586 , BY265614 , BY293803
UniProt ID:	Q9DB32
Synonyms:	1500017E18Rik; 2810014I23Rik; C330022E15Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)