

Product datasheet for **MP206395**

Hsd17b4 Mouse qPCR Primer Pair (NM_008292)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15488
Forward Sequence:	TGGCTTTGCCATGAGAGCTGTG
Reverse Sequence:	CTGATTCCGCTTTCTGACGATGG
ACCN:	BC022175 , NM_008292 , NM_008292.1 , NM_008292.2 , NM_008292.3 , NM_008292.4
UniProt ID:	P51660
Synonyms:	17-beta-HSD; 17[b]-HSD; DBP; MFE-2; Mfp-2; MFP2; MPF-2; perMFE-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)