

## Product datasheet for **MP206461**

### Icosl Mouse qPCR Primer Pair (NM\_015790)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50723
Forward Sequence:	TCTTGAAGAGGTGGTCAGGCT
Reverse Sequence:	GCCATTCTGGACATGCAGGTG
ACCN:	<u><a href="#">BC029227</a></u> , <u><a href="#">NM_015790</a></u> , <u><a href="#">NM_015790.1</a></u> , <u><a href="#">NM_015790.2</a></u> , <u><a href="#">NM_015790.3</a></u>
UniProt ID:	<u><a href="#">Q9JHJ8</a></u>
Synonyms:	AU044799; B7-H2; B7h; B7RP-1; BG071784; GI50; GL50; GL50-B; ICOS-L; Icoslg; LICOS; Ly115I
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)