

## Product datasheet for **MP206467**

### **Ikzf2 Mouse qPCR Primer Pair (NM\_011770)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	22779
<b>Forward Sequence:</b>	CCAGAATGTCAGCATGGAGGCT
<b>Reverse Sequence:</b>	AGCAGGTCTCTCAAAGGCACC
<b>ACCN:</b>	<u><a href="#">NM_011770</a></u> , <u><a href="#">NM_011770.1</a></u> , <u><a href="#">NM_011770.2</a></u> , <u><a href="#">NM_011770.3</a></u> , <u><a href="#">NM_011770.4</a></u> , <u><a href="#">BC138606</a></u> , <u><a href="#">BC050049</a></u> , <u><a href="#">BC055378</a></u> , <u><a href="#">BC138608</a></u>
<b>UniProt ID:</b>	<u><a href="#">P81183</a></u>
<b>Synonyms:</b>	A730095J18Rik; Helios; Zfpn1a2; Znfn1a2
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)