

Product datasheet for **MP206507**

Ipo9 Mouse qPCR Primer Pair (NM_153774)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226432
Forward Sequence:	CAGGCAACAGTCAGTGGTCTTC
Reverse Sequence:	GCTGAAGCACATGCGTACTCTC
ACCN:	<u>BC098508</u> , <u>NM_153774</u> , <u>NM_153774.1</u> , <u>BC048549</u> , <u>BC056469</u> , <u>BC079546</u>
Synonyms:	0710008K06Rik; AI845704; C78347; Imp9; mKIAA1192; Ranbp9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)