

Product datasheet for **MP206514**

Icos Mouse qPCR Primer Pair (NM_017480)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54167
Forward Sequence:	GCAGCTTTCGTTGTGGTACTCC
Reverse Sequence:	TGTGTTGACTGCCCCATGAAC
ACCN:	BC034852 , NM_017480 , NM_017480.1 , NM_017480.2 , BB850356
UniProt ID:	Q9WVS0
Synonyms:	AILIM; CCLP; CRP-1; H4; Ly115
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)