

Product datasheet for **MP206526**

Ing1 Mouse qPCR Primer Pair (NM_011919)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26356
Forward Sequence:	TGCAGAGGAACGTCTCGCTGAT
Reverse Sequence:	TGCCGTCTGTCTCCGTTTGAA
ACCN:	<u>NM_011919</u> , <u>NM_011919.1</u> , <u>NM_011919.2</u> , <u>NM_011919.3</u> , <u>NM_011919.4</u> , <u>NM_011919.5</u> , <u>BC147770</u> , <u>BC016573</u> , <u>BC147784</u> , <u>BC152847</u> , <u>BX528243</u> , <u>BY063246</u>
UniProt ID:	<u>Q9QXV3</u>
Synonyms:	2610028J21Rik; AA407184; A1875420; mING1h; p33Ing1; p37Ing1b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

