

## Product datasheet for **MP206554**

### **Bpnt2 Mouse qPCR Primer Pair (NM\_177730)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	242291
<b>Forward Sequence:</b>	CCGCAAGATGTTCTACCTGCTC
<b>Reverse Sequence:</b>	GGATGTCCTCAGGAATCTCCG
<b>ACCN:</b>	<u><a href="#">BC048776</a></u> , <u><a href="#">NM_177730</a></u> , <u><a href="#">NM_177730.1</a></u> , <u><a href="#">NM_177730.2</a></u> , <u><a href="#">NM_177730.3</a></u> , <u><a href="#">BC138209</a></u> , <u><a href="#">BC145952</a></u> , <u><a href="#">NM_177730.4</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q80V26</a></u>
<b>Synonyms:</b>	1110001C20Rik; AA408880; AI451589; AL022796; B230207P20; gPAPP; Jaws
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)