

## Product datasheet for **MP206566**

### Ide Mouse qPCR Primer Pair (NM\_031156)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	15925
Forward Sequence:	CAAACCTCTCCTTCCAAGTCAGC
Reverse Sequence:	TGTTCTCCGAGGTGCTCTGCAT
ACCN:	<u><a href="#">BC041675</a></u> , <u><a href="#">NM_031156</a></u> , <u><a href="#">NM_031156.1</a></u> , <u><a href="#">NM_031156.2</a></u> , <u><a href="#">NM_031156.3</a></u> , <u><a href="#">NM_031156.4</a></u>
Synonyms:	1300012G03Rik; 4833415K22Rik; AA675336; AI507533
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)