

Product datasheet for **MP206616**

Itga3 Mouse qPCR Primer Pair (NM_013565)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16400
Forward Sequence:	ACGAGACCATCCTGTGTGAGCT
Reverse Sequence:	GCAAGTTGTCCTGGTGACTCGA
ACCN:	<u>BC053031</u> , <u>BC062205</u> , <u>NM_013565</u> , <u>NM_013565.1</u> , <u>NM_013565.2</u> , <u>NM_013565.3</u> , <u>BC023259</u> , <u>BC096741</u> , <u>BQ032311</u>
UniProt ID:	<u>Q62470</u>
Synonyms:	AA407068; CD49C; GAPB3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)