

Product datasheet for **MP206627**

Itgam Mouse qPCR Primer Pair (NM_001082960)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16409
Forward Sequence:	TACTTCGGGCAGTCTCTGAGTG
Reverse Sequence:	ATGGTTGCCTCCAGTCTCAGCA
ACCN:	NM_001082960 , NM_001082960.1 , BC126892 , BC156094 , BC156991
Synonyms:	Cd11b; CD11b/CD18; CR3; CR3A; F730045J24Rik; Ly-40; Mac-1; Mac-1a; MAC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)