

Product datasheet for **MP206640**

Itgb6 Mouse qPCR Primer Pair (NM_021359)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16420
Forward Sequence:	ACTCATTCTGGAGCAACCGTG
Reverse Sequence:	GCTGTGAAAGACAGGTTGAGTCC
ACCN:	<u>BC049185, NM_021359, NM_021359.1, NM_021359.2, NM_021359.3</u>
UniProt ID:	<u>Q9Z0T9</u>
Synonyms:	2210409C20Rik; 4831415H04Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)