

Product datasheet for **MP206650**

Ilkap Mouse qPCR Primer Pair (NM_023343)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67444
Forward Sequence:	TCTCTTGCCACATCAGGTTCCC
Reverse Sequence:	CTGCCACATAGCCTTTCAGACC
ACCN:	<u>BC026953</u> , <u>BC027439</u> , <u>NM_023343</u> , <u>NM_023343.1</u> , <u>NM_023343.2</u> , <u>BC027439.1</u> , <u>BC011427</u>
UniProt ID:	<u>Q8R0F6</u>
Synonyms:	0710007A14Rik; 1600009O09Rik; PP2C-DELTA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)