

## Product datasheet for **MP206749**

### Il13ra1 Mouse qPCR Primer Pair (NM\_133990)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16164
Forward Sequence:	TGTTCTTGCCGACGCTGTCTAC
Reverse Sequence:	GTGTAGAAGGTGGAGTTTTGCTC
ACCN:	<u><a href="#">BC052425</a></u> , <u><a href="#">BC059939</a></u> , <u><a href="#">NM_133990</a></u> , <u><a href="#">NM_133990.1</a></u> , <u><a href="#">NM_133990.2</a></u> , <u><a href="#">NM_133990.3</a></u> , <u><a href="#">NM_133990.4</a></u> , <u><a href="#">NM_133990.5</a></u> , <u><a href="#">BC021472</a></u>
UniProt ID:	<u><a href="#">O09030</a></u>
Synonyms:	AI882074; CD213a; CD213a1; IL-13R-alpha-1; IL-13r[a]; Il13ra; NR4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)