

Product datasheet for **MP206750**

IL13ra2 Mouse qPCR Primer Pair (NM_008356)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16165
Forward Sequence:	TGGTCTGCTCTTGAAACCTGG
Reverse Sequence:	GCTGGAGGTAATCAGCACACTG
ACCN:	<u>BC003723</u> , <u>NM_008356</u> , <u>NM_008356.1</u> , <u>NM_008356.2</u> , <u>NM_008356.3</u> , <u>NM_008356.4</u>
UniProt ID:	<u>O88786</u>
Synonyms:	CD213a; CD213a2; IL-13R-alpha-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)