

Product datasheet for **MP206756**

Il17rc Mouse qPCR Primer Pair (NM_134159)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171095
Forward Sequence:	TAGAGCCAGACTCTGAGAGGGT
Reverse Sequence:	AAGGCGCATCTAGCTGCCATAC
ACCN:	<u>BC004759</u> , <u>NM_134159</u> , <u>NM_134159.1</u> , <u>NM_134159.2</u> , <u>NM_134159.3</u> , <u>NM_134159.4</u> , <u>BQ551040</u> , <u>BY018830</u>
UniProt ID:	<u>Q8K4C2</u>
Synonyms:	1110025H02Rik; Gm19850; IL-17RC; IL17-RC; IL17-RL; Il17rl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)