

Product datasheet for **MP206771**

Il2rb Mouse qPCR Primer Pair (NM_008368)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 16185 |
| Forward Sequence: | CTCAAGTGCCACATCCCAGATC |
| Reverse Sequence: | AGCACTTCCAGCGGAGAGATCT |
| ACCN: | BC026869 , BC099616 , NM_008368 , NM_008368.1 , NM_008368.2 , NM_008368.3 , NM_008368.4 , BC028446 , BQ033566 |
| UniProt ID: | P16297 |
| Synonyms: | CD122; Il-2/15Rbeta; Il-2R; Il-2Rbeta; IL-15R; IL-15Rbeta; IL15R; IL15Rbeta; p70 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)