

Product datasheet for **MP206810**

Irak1bp1 Mouse qPCR Primer Pair (NM_022986)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	65099
Forward Sequence:	CTCTGTTGTCATCAGCCACCT
Reverse Sequence:	CAACGAGGTCACAGACTTCCTG
ACCN:	<u>BC061099</u> , <u>NM_022986</u> , <u>NM_022986.1</u> , <u>NM_022986.2</u> , <u>NM_022986.3</u> , <u>NM_022986.4</u> , <u>BC046301</u> , <u>BY719456</u>
UniProt ID:	<u>Q9ESJ7</u>
Synonyms:	4921528N06Rik; Aabp3; Al851240; Aip70; Simpl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

