

## Product datasheet for **MP206813**

### Irak4 Mouse qPCR Primer Pair (NM\_029926)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	266632
Forward Sequence:	GGTCCAGATTGAGCTGTTTGCC
Reverse Sequence:	GTTTGTGCCACTGTTGCCGCTT
ACCN:	<a href="#">BC051676</a> , <a href="#">NM_029926</a> , <a href="#">NM_029926.1</a> , <a href="#">NM_029926.2</a> , <a href="#">NM_029926.3</a> , <a href="#">NM_029926.4</a> , <a href="#">NM_029926.5</a> , <a href="#">BE635313</a>
UniProt ID:	<a href="#">Q8R4K2</a>
Synonyms:	8430405M07Rik; 9330209D03Rik; IRAK-4; NY-REN-64
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)