

Product datasheet for **MP206819**

Itsn2 Mouse qPCR Primer Pair (NM_011365)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20403
Forward Sequence:	CCAGTTACCTCCACAGCCTATC
Reverse Sequence:	CCACTCTCCATCTTTCTGGGTC
ACCN:	<u>NM_011365</u> , <u>NM_011365.1</u> , <u>NM_011365.2</u> , <u>NM_011365.3</u> , <u>NM_011365.4</u> , <u>BC037119</u> , <u>BC138262</u> , <u>BC171950</u> , <u>BE650375</u> , <u>BI872589</u>
Synonyms:	AI327390; Ese2; mKIAA1256; Sh3d1B; Sh3p18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)