

## Product datasheet for **MP206835**

### lyd Mouse qPCR Primer Pair (NM\_027391)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70337
Forward Sequence:	CAACGAGATCAGTGTGCCATCG
Reverse Sequence:	TCAGACGAGGACCACAGTTGAG
ACCN:	<u><a href="#">BC023358</a>, <a href="#">NM_027391</a>, <a href="#">NM_027391.1</a>, <a href="#">NM_027391.2</a>, <a href="#">NM_027391.3</a>, <a href="#">NM_027391.4</a></u>
UniProt ID:	<u><a href="#">Q9DCX8</a></u>
Synonyms:	0610009A07Rik; AI265638; Dehal1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)